

Microsoft Excel® Data Acquisition Module

OVERVIEW-

Abstract

The PLA 3.0 Data Acquisition Module for Microsoft Excel enables users to import data from Excel sheets into biological assays for further analysis.

Features

- Fast, secure and reliable data import from Excel files.
- Option to save the configuration to repeat the process with the same settings.
- All imported data is fully retraceable (audit trail).
- Fully customizable for import from different sheets, regions.

More information on PLA 3.0 Data Acquisition Modules:

TECHNICAL SPECIFICATION

Type

Data Acquisition Module

Version

1.2.0 (Build 99)

System Requirements

- PLA 3.0.4
- Bioassay Package 23 or higher

Supported Systems

Microsoft Excel 97 and newer

Supported File Formats

- xls
- .xlsx

Supported Languages

English, German

DELIVERY AND SUPPORT-

Release Date

2016/03/25

Product No.

24032

Scope of Delivery

- Software (CD or download)
- User manual (PDF)

Lifecycle

We will provide support for PLA 3.0 and all related modules until March 15, 2024.

Support and Upgrades

This module is covered by our PLA Support Contract. The first year is included and the contract can be renewed year by year. To find out more about its benefits, please go to:




SCREENSHOTS -


Creating a data acquisition configuration.



Managing the settings.



Editing the configuration.



Successful acquisition of Excel data.

